





Consumable price list




***B = Box**



***S = Single**



***P = Packet**



Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
500mL Sterile Dosed (MICRO) 	TVC (22 & 3C) Total coliforms E coli Enterococci Pseudomonas spp. Or Pseudomonas aeruginosa Sulphite Reducing Bacteria Legionella	Supplied by manufacturer with sodium thiosulphate added. Must be used before expiry date shown on bottle. Do not use if seal is broken	Do not rinse Fill to shoulder Do not overflow	Yes Legionella – No Other microbes –Yes	Micro	£56.32 (144pcs)	£0.39	Consum1
1Lt Sterile Dosed (MICRO) 	TVC (22 & 3C) Total coliforms E coli Enterococci Pseudomonas spp. Or Pseudomonas aeruginosa Sulphite Reducing Bacteria Legionella	Supplied by manufacturer with sodium thiosulphate added. Must be used before expiry date shown on bottle. Do not use if seal is broken	Do not rinse Fill to shoulder Do not overflow	Yes Legionella – No Other microbes –Yes	Micro	£40.95	£0.59	Consum2



Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
1LT Not Sterile, Undosed – (Chemistry) 	Chemistry tests Metals Anions/cations pH Color turbidity Effluent Nitrogen tests	none	Do not rinse Fill to shoulder	none	Chemistry	£20.95 (70pcs)	£0.30	Consum3
10L Grab Disposable 	Cryptosporidium Giardia	None	Fill to shoulder Rinse with sample before filling	No	Subcontract	Available upon request	£5	Consum4



Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
<p>Blue Sponge with 10ml Neutralising buffer</p> 	Good for presence / absence tests, not suitable for quantification tests	Neutralizing buffer	Open outer packaging, take swab aseptically using the bag to hold the swab. Return, in sealed inside bag.	Yes	Micro	£162.00 (200pcs)	£0.81	Consum5
<p>Pathogen Swab</p> 	Good for presence / absence tests, not suitable for quantification tests e.g., Listeria, Salmonella	Neutralizing buffer	Remove swab, swab target area, replace swab into tube and secure.	Yes	Micro	£116.00 (100pcs)	£2.00	Consum6
<p>EnviroStik – blue sponge with straight handle with 10ml Neutralising Buffer</p> 	<p>Good for presence / absence tests, not suitable for quantification tests</p> <p>Listeria swabs</p>	Neutralizing buffer	Open outer packaging, take swab aseptically by only holding the handle. Return, in sealed inside bag.	Yes	Micro	£154.49 (125pcs)	£1.24	Consum7


Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
TS/5-42 10ML Neutralising buffer 	Good for enumeration tests (e.g., TVC, enteros, Staph) can be used for presence / absence tests. Several enumerations tests can be requested from one swab, but only one presence / absence test per swab.	Neutralizing buffer	Remove swab, swab target area, replace swab into tube and secure.	Yes	Micro	£34.74 (100pcs)	£0.35	Consum8
150ML Red cap 	Sterile 150ml container for food and liquids. Does not contain sodium thiosulphate so not suitable for most water testing. Taste and Odour Chemistry	None	Fill to line Do not overflow	Depends on matrix, frozen, chilled or ambient	Micro & Chem	£46.04 (350pk)	£0.18	Consum9

Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
Dry Gluten, Soya & Egg Swab – TS/17-A-SLS – IN HOUSE only Testing 	Gluten, Soya & Egg	None	None	No	Chemistry	£39.74 (250pcs)	£0.16	Consum10
S02 Swab- in house testing 	Sulphur dioxide	None	Remove swab, swab area to be tested return to tube and close lid. Do not discard contents	No	Chemistry	£128.52 (50pcs)	£2.57	Consum11

Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
DNA 	Speciation tests Chicken Turkey Donkey Horse Goat Sheep Pig Bovine/Beef Haddock Cod	None	No	No	Subcontract	Available upon request	£2 per swab	Consum12
Bio Check swab 	Milk, Egg, Soya, Protein – not for inhouse testing	None	No	No	Subcontract	Available upon request	£2 per swab	Consum13

Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
300ml Amber Glass/1litre glass 	PAHs Cyanides Lead (1 litre is available for fat/oil and grease testing)	None	Do not rinse Fill to shoulder	No	Subcontract	Available upon request	Available upon request	Consum14
Blender Bag 	n/a	n/a	n/a	n/a	n/a	n/a	£51.60 (box of 500)	Consum15

Description	Parameters Analyzed	Preservation	Filling Instructions	Refrigerate?	Area of use	Price per box	Individual Cost	ID Code
Temperature Bags (RSA clients) They are now white 	For clients that supply products to own brand Tesco and M&S	n/a	n/a	n/a	n/a	n/a	n/a	Consum16
Grip Bags 	8"X11"	n/a	n/a	n/a	n/a	£83.77 Pk1000	£8.38 per 100 pack	Consum17

<p>BARCODE LABELS</p> 	n/a	n/a	n/a	n/a	n/a	N/A	N/A	Consum18
<p>AIRPLATES</p>	<p>We will need to know. No. of TVC: _____ No. of Y&M: _____</p>					N/A	£0.77	Consum19